
FIBRAPOST PLUS

Glass Fibre reinforced Composite Posts for the corono-radicular dental reconstruction.

Directions for use

This instruction manual refers to the use of SEALACORE DC (PD; Switzerland) products. **Note:** for the use of another bonding and/or core build up system read the instructions of the manufacturer.

Endodontically prepare the canal in the usual manner, clean and dry. Obtainate the apical portion.

- 1 Choose the post size and measure the length of the canal.
- 2-3 Drill the canal with the corresponding Drill-A-Post (800-1'200 RPM), placing its depth stopper at the required length.
- 4 Check the size of the post and its position in the canal.
- 5 Excess in length may be adjusted with a diamond disc. **Carry it out outside of the mouth**, under aspiration. Clean the post with alcohol. No need to silanize (FIBRAPOST PLUS are pre-silanized).
- 6 SEALACORE DC Nano-reinforced Self-Etch Bond. Mix 1 drop of Liquid A and 1 drop of Liquid B on a mixing pad with the application brush for approximately 2 seconds. Avoid exposure to intensive light. Apply a layer with the application brush to the enamel/dentine and rub into the tooth surface for 20 seconds.
- 7 Dry the adhesive layer for at least 5 seconds with an air syringe.
- 8 Do not light cure at this step of the procedure.
- 9 Apply SEALACORE DC Composite directly into the root canal using thin intraoral tips. **Note: Do not light cure at that stage!**
- 10 Insert the FIBRAPOST making sure to leave a slight coronal excess for further anchorage of the core.
- 11 Light-cure for at least 40 seconds.
- 12 **Core build up.** Place a matrix around the prepared tooth and apply SEALACORE DC Composite directly from the cartridge/double-barrel syringe. Place the tip of the polymerization lamp as close as possible to the surface of the composite and polymerize with halogen light for 40 seconds. In case of a larger core build-up, generally apply in layers of 2 mm and cure each layer for 10 seconds.
Final Polymerisation: Light cure 40 seconds or wait 5 min. After curing, the composite can be worked immediately.

Properties

FIBRAPOST PLUS are made from specially treated Glass Fibres (approx. 65%) reinforced with, and firmly bound to, an UDMA (Urethane Dimethacrylate) resin (approx. 20%) known for its resistance, safety and capacity to bind to Composite resins used for post cementation.

Features and benefits

- ▶ Optimal flexural strength (~ 35 GPa) close to the one of the dentine
- ▶ High resistance (~ 1'300 MPa) thanks to interlaminar fibres cohesion
- ▶ Natural shade for aesthetic reconstruction
- ▶ Convenient Radioopacity (~ 200% Al (Ø 2.5 mm))
- ▶ High adhesion capacity between post and cement
- ▶ Biocompatibility

Indication

Aesthetic post-retained restorations and core building-up of endodontically treated non-vital teeth.

Sizing and shape

4 sizes according to canal size allowing the preservation of maximum healthy dentine. Cylindro-conical shape.

Contraindications

Do not use when root is not firmly anchored in the bone, or in case of insufficient tooth structure available. Allergy to methacrylates.

Disinfection

Disinfect with alcohol before using. Single use.

Storage Conditions

5°C – 30°C. Protected from humidity.

Packaging

ALL-IN-ONE BUILD UP KIT
FIBRAPOST PLUS
4 x 4 glass fiber composite posts (sizes 1, 2, 3 & 4)
4 drills with adjustable depth stopper
SEALACORE DC
1 dual-cure nano-reinforced self-etch bond (liquid A 4ml + liquid B 4 ml)
1 dual-cure composite (syringe 5 ml)
1 mixers
10 tips
25 minibrushes
1 mixing plate

REF 10475

Assortment classic

6 posts each N° 1 to 4 = 24 posts
1 reamer each N° 1 to 4 = 4 reamers

REF 22595

Refills (in racks of 12 pieces)

N° 1 22596
N° 2 22597
N° 3 22598
N° 4 22599
N° 1 to 4 (3 pieces of each size) 22600

Size	N°1	N°2	N°3	N°4
(L) Total length in mm	19	19	19	19
(l) Length of the tapered part in mm	11	11	11	11
(D1) Ø in mm at the tip	0.85	1.06	1.26	1.47
(D2) Ø in mm at 11 mm	1.30	1.50	1.70	1.90

For dental professional use. Made in Switzerland.

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